



OrthoSleeve[®] **Braces and Supports**

World's thinnest,
most comfortable,
easy-to-wear &
effective braces



WORK

PLAY

LIVE





making life better.



The Difference...

Compression Zone Technology® – It's What's Up Our Sleeves

A lot of our delighted customers want to know what it is that makes OrthoSleeve bracing sleeves so different from all the compression stockings, splints, and braces used to treat painful conditions like Plantar Fasciitis, Patellar tendonitis, knee pain, tennis elbow and Carpal Tunnel Syndrome.

Traditional splints and braces; Patient compliance challenges

Traditional treatments like splints and braces rely on the patient to tighten velcro straps or other similar features to provide adequate support to the underlying tissue. While they certainly have some benefit, there's a downside: the patient has no way of knowing if the way they are tightening the product is affecting their circulation until it's already been on for a long time and they notice symptoms like discoloration or temperature changes in the skin around the brace.

It's not easy to strike that perfect balance in which the brace is tight enough to function properly, but not too tight to harm circulation.

In addition, most of these braces and splints are made of materials that are hot, heavy, bulky, and generally uncomfortable for long-term use (such as overnight or throughout a work day).

The end result is that the device isn't worn as often as it should be to treat the condition from which the patient is suffering.

The CZT® difference

If you've tried on an OrthoSleeve product, you noticed that each sleeve offers medical grade orthopedic support like a strap or brace, but unlike these other products, an OrthoSleeve is thin, light, comfortable, and made of moisture-wicking fabric that's easy to slip on and off.

So what makes the difference?

The design of OrthoSleeve bracing sleeves is based on a patented formula called Compression Zone Technology®.

Without going too deep into the science behind the technology, picture each OrthoSleeve bracing sleeve divided into a number of zones with each zone providing a different level of compression for orthopedic effect. Generally, the lightest compression is at the outer edges and openings of the sleeve, which facilitates the easy-on-easy-off design and promotes free circulation of blood into and out of the covered area. Then, as you move down the sleeve from either end, the zones of compression change to accommodate the part of the body being treated, with the firmest compression to be found at the exact spot where the highest level of orthopedic support is required to treat the condition in question.

All of this is accomplished through the use of high-grade compression fibers woven together seamlessly to provide maximum benefit and maximum comfort in the same thin and light product, without resorting to any straps or patient-controlled tension.

Why Compression Zone Technology® is better

The bottom line is comfort. By creating a product that is highly effective in providing pain relief, promoting circulation, and supporting injured tendons or muscles, but that is comfortable enough to wear all day or night even under your normal clothing, OrthoSleeve has offered millions of sufferers the best option available:

Comfortable, mobile pain relief you can wear any time of day or night, so you can stay on the go.

WWW.ORTHOSLEEVE.COM

OrthoSleeve



The Orthosleeve family of products combines the benefits of medical grade orthopedic support with the power of graduated compression in one thin, light, comfortable sleeve that relieves pain and improves mobility for millions of people. Compression Zone Technology® built into every OrthoSleeve product is scientifically precise graduated compression that surrounds your joints, muscles, tendons, and ligaments in anatomically correct zones to improve circulation, support and stabilize joint structure, and reduce swelling. The result? Pain relief, increased mobility, greater confidence and comfort as you pursue your active lifestyle.

OrthoSleeve Compression Zone Technology® (CZT®) defines a new way of making ortho supports. This breakthrough technology delivers medical grade support in thinner, more comfortable, discreet medical devices. No longer is it necessary to use heavy splints, uncomfortable neoprene fabrics and frustrating Velcro.

Introducing K-Zone® Technology (KZT®) a breakthrough design that provides the support of a kinesiology taping treatment without the inconvenience of its temporary and binding nature. Using a hypoallergenic gel stabilizer system with high resistance fabrication zones, KZT creates a strapping type effect in specific areas to steady and help secure supportive joint structures. The result is stabilization without immobilization. Proper medical-grade compression in combination with KZT provides a slip-on therapy solution for use whenever you need it. KZT is found in the newest OrthoSleeve product, the AF7 Ankle Bracing Sleeve.

Orthosleeve products allow you to receive the support and pain relief you deserve any time of day. You can sleep in your OrthoSleeves to awake pain-free. You can easily slip on your OrthoSleeve under any type garment... business suit, work clothes, gym clothes, socks...you name it. These products go where you go. We call this *motion comfort™*.

This Compression Zone Technology® uses patented graduated compression for vascular improvement or vascular safety. OrthoSleeve products give you circulation improvement and support without circulation interference, unlike most bulky orthopedic and sports supports.





RELIEVES & PREVENTS:

- PLANTAR FASCIITIS
- ACHILLES TENDONITIS
- HEEL PAIN

AIDS RECOVERY
REDUCES SWELLING

The Patented FS6 Foot Bracing Sleeve has sold millions of units worldwide because it's that incredibly effective at relieving the pain of Plantar Fasciitis, chronic arch and heel pain, and swelling in the foot. Compression Zone Technology® woven into this amazing foot sleeve provides medical grade orthopedic support for the foot structure and the Achilles tendon, graduated compression to improve circulation and reduce swelling, all in an ultra-comfortable moisture-wicking foot sleeve you can wear under your normal socks and shoes. Wear FS6 overnight and experience morning pain relief from the first day.



TARGETED CONDITIONS

- Plantar Fasciitis
Achilles Tendonitis
Swollen Feet
Arch Pain
Heel Pain
Sever's Disease



- ZONE 6**
Light Compression
For a non-binding fit
- ZONE 5**
Moderate Compression
To support the Achilles, lift the Plantar Fascia & boost venous system
- ZONE 4**
Moderate Compression
With Shaped-To-Fit Design
- ZONE 3**
Firm Compression
To support venous-flow & lift the Plantar Fascia
- ZONE 2**
Moderate Compression
To boost venous-flow & support foot structure
- ZONE 1**
Light Compression
Wide, smooth, flat opening

MILLIONS SOLD!

Tested and Proven Effective
Excerpt from the FS6 white paper by Roy H. Lidtke DPM, CPED, FACFAOM
“The FS6 Day or Night Compression Foot Sleeve provides graduated compression across the arch, foot and ankle at a level to reduce edema and the pain associated with conditions such as Plantar Fasciitis, heel spurs, and other related conditions. The dynamic gait adaptation seen with the FS6 sleeves suggests that it may be producing a neuromuscular adaptation to improve foot function.”



Where to measure	Average Shoe Size		
	Size	Arch Circumference	Average Shoe Size
	XS	4-6 in/8-16 cm	US: Y12-M3/Y12-W4
	S	5-8 in/13-20 cm	US: M 3-5.5/W 4-6.5 UK: 1-3 EU: 33-36
	M	7-10 in/18-25 cm	US: M 6-9/W 7-10.5 UK: 4-8 EU: 37-41
	L	9-13 in/23-33 cm	US: M 10-13/W 11-Plus UK: 9-13 EU: 42-46
	XL	11-15 in/28-38 cm	US: M 13.5- 15 UK: 14-16 EU: 47-50
	XXL	13-18 in/33-46 cm	US: M 15.5-PLUS UK: 17-Plus EU 51-PLUS

	32339B	XS		32341W	S
	32341B	S		32343W	M
	32343B	M		32345W	L
	32345B	L		32347W	XL
	32347B	XL			
	32349B	XXL			
	32339N	XS		32341PF	S
	32341N	S		32343PF	M
	32343N	M		32345PF	L
	32345N	L		32347PF	XL
	32347N	XL			
				32341RY	S
				32343RY	M
				32345RY	L
				32347RY	XL





PROVIDES MEDIAL &
LATERAL STABILIZATION
AIDS RECOVERY
REDUCES SWELLING
IMPROVED BALANCE
FALL PREVENTION

The Patent Pending AF7 Ankle Bracing Sleeve is a one-of-a-kind ankle brace designed to stabilize without immobilizing, for relief and balance while active. It provides medial and lateral stabilization by combining Compression Zone Technology® with our patent pending K-Zone® Technology (KZT®). KZT keeps the bracing sleeve in position while providing extra support and stabilization where needed most. With a thin profile, the AF7 fits under any sock and inside any shoe making it more comfortable and effective than traditional splints and braces. Designed to fit individual shape and size.

TARGETED CONDITIONS

Ankle Pain
Weak Ankles
Achilles Tendonitis



- ZONE 7**
Light Compression
Wide, smooth band for comfort & ease of donning
- ZONE 6**
Moderate Compression
For increased circulation & ankle stability
- ZONE 5**
Firm Compression
For increased circulation & ankle joint control
- ZONE 4**
Moderate Compression
Anatomical "Pressure Relief Zone"
- ZONE 3**
Firm Compression
For increased circulation & ankle joint control
- ZONE 2**
Moderate Compression
Boost venous flow & stabilize
- ZONE 1**
Light Compression
Wide, smooth, flat opening

Introducing K-Zone® Technology (KZT®)

The KZT® Ankle Stabilizer Zone acts as a woven "stay" developed using high resistance fabrication techniques. This eliminates the use of a plastic or metal bar for stabilizing the ankle and, instead, creates a flexible rigidity for support. This Ankle Stabilizer Zone fuses into the Lateral Gel Stabilizer Bar for a seamless bracing effect similar to various clinical taping techniques.

The KZT® Lateral Gel Stabilizer Bar is a one-of-a-kind bracing element designed to stabilize the Talus and midtarsal joints without immobilizing. This provides lateral stabilization for inversion resistance to neutralize unsafe motion. By combining KZT® with Compression Zone Technology® optimal positioning is achieved, for multiple type ankle injury prevention.



Where to Measure	Size	Ankle Circumference	Average Shoe Size
	S	6.5-8.5 in/17-22 cm	US: M 5-7/W 6-8 UK: 2-4 EU: 35-38
	M	8.5-10.5 in/22-27 cm	US: M 8-9/W 9-10 UK: 5-9 EU: 39-42
	L	9.5-12.5 in/24-32 cm	US: M 10-12/W 11-13 UK: 10-12 EU: 43-46
	XL	12-14 in/30-36 cm	US: M 13-15 UK: 13-17 EU: 47-50
	XXL	13-16 in/33-40 cm	US: M 15-Plus UK: 17-Plus EU: 51-Plus

Black-Left Foot

35340B	S
35342B	M
35344B	L
35346B	XL
35348B	XXL

Black-Right Foot

35350B	S
35352B	M
35354B	L
35356B	XL
35358B	XXL





IMPROVES CIRCULATION
EASES MUSCLE SORENESS
RELIEVES TIRED, ACHY LEGS
AIDS RECOVERY
EASY DONNING

The Patented CS6 Calf Bracing Sleeve offers a comfortable and easy-to-wear alternative to tight, binding compression stockings. Compression Zone Technology® uses specific graduated compression to encourage efficient blood flow and oxygen return, providing relief from common conditions like chronic venous insufficiency and edema. It also provides medical grade support for the Achilles tendon and calf muscles, relieving shin splints, cramping, and Achilles tendonitis. Fits easily under your normal socks and clothing and can replace compression stockings.

TARGETED CONDITIONS

Venous Insufficiency
Poor Circulation
Varicose Veins
Achilles Tendonitis
Restless Leg Syndrome
Aching Legs & Feet
Leg Cramps
Shin Splints



ZONE 6
Light Compression
Gentle compression, wide smooth hem for non-binding comfort

ZONE 5
Moderate Compression
Graduated compression with "Calf Stabilizer Zone" for injury protection

ZONE 4
Moderate Compression
Graduated compression for muscle support & continued circulation support

ZONE 3
Firm Compression
Firm graduated compression for ligament support & improved circulation

ZONE 2
Firm Compression
Firm Medical Grade compression for circulation boost & faster recovery with "Achilles Stabilizer Zone"

ZONE 1
Light Compression
Gentle compression with smooth flat opening that comfortably fits under socks

Where to measure	Size	Ankle Circumference	Calf Circumference	Average Shoe Size
	S	6.5-8.5 in/17-22 cm	9.5-14.5 in/25-37 cm	US: M 3-5.5/W 4-6.5 UK: 1-3 EU 33-36
	M	8.5-10.5 in/22-27 cm	12-17 in/30-43 cm	US: M 6-9.5/W 7-10.5 UK: 4-8 EU: 37-41
	L	9.5-12.5 in/24-32 cm	14-19 in/35-48 cm	US: M 10-13/W 11-Plus UK: 14-16 EU: 43-47
	XL	12-14 in/30-36 cm	14.5-20 in/37-50 cm	US: M 13.5-15 UK: 14-16 EU: 47-50

Black	42341B 42343B 42345B 42347B	S M L XL	Pink Fusion	42341PF 42343PF 42345PF	S M L
White	42341W 42343W 42345W 42347W	S M L XL	Natural	42341N 42343N 42345N	S M L
			Reflector Yellow	42341RY 42343RY 42345RY 42347RY	S M L XL

FS6+ Foot + Calf Bracing Sleeve

RELIEVES PLANTAR FASCIITIS
EASES HEEL PAIN
REDUCES SWELLING
IMPROVES CIRCULATION
EASES MUSCLE SORENESS
RELIEVES TIRED, ACHY LEGS
AIDS RECOVERY

The Patented FS6+ Foot + Calf Bracing Sleeve provides medical grade orthopedic support and graduated compression to ease the pain of Plantar Fasciitis, Achilles tendonitis, and more using Compression Zone Technology®. The FS6+ is a new innovative design which combines the benefits of the FS6 and CS6 bracing sleeves in one easy to use orthopedic support. Thin, light and more comfortable to wear, the FS6+ offers a unique, premium alternative to traditional compression socks & hosiery that can make you sweat. Our moisture-wicking fabric keeps your legs and feet cool throughout the day.

TARGETED CONDITIONS

Plantar Fasciitis
Achilles Tendonitis
Aching Legs & Feet
Swollen Feet
Arch Pain
Heel Pain
Leg Cramps
Shin Splints

ZONE 11

Light Compression
With smooth, flat closure

ZONE 9 & 10

Light Compression
Graduated compression with "Calf Stabilizer Zone" for injury protection

ZONE 8

Moderate Compression
Graduated compression for muscle support & continued circulation support

ZONE 7

Firm Compression
Graduated compression for ligament support & improved circulation

ZONE 5 & 6

Firm Compression
To support the Achilles, lift the Plantar Fascia

ZONE 4

Moderate Compression
With shaped-to-fit comfort

ZONE 3

Firm Compression
To support venous-flow and lift the Plantar Fascia

ZONE 2

Moderate Compression
To boost venous-flow and support foot structure

ZONE 1

Light Compression
With smooth, flat closure



Where to measure	Size	Arch Circumference	Ankle Circumference	Calf Circumference	Average Shoe Size
	S	5-8 in/13-20 cm	6.5-8.5 in/17-22 cm	9.5-14.5 in/24-37 cm	US: M3-5.5/W 4-6.5 UK: 1-3 EU: 33-36
	M	7-10 in/18-25 cm	8.5-10.5 in/22-27 cm	12-17 in/30-43 cm	US: M 6-9.5/W 7-10.5 UK: 4-8 EU: 37-41
	L	9-13 in/23-33 cm	9.5-12.5 in/24-32 cm	14-19 in/35-48 cm	US: M 10-13/W 11-Plus UK: 9-13 EU: 42-46
	XL	11-15 in/28-38 cm	12-14 in/30-36 cm	14.5-20 in/37-51 cm	US: M 13.5-15 UK: 14-16 EU: 47-50
	XXL	13-18 in/33-46 cm	13-16 in/33-40 cm	16-23 in/40-58 cm	US: M 15.5-Plus UK: 17-Plus EU 51-Plus

Black

22341B
22343B
22345B
22347B
22349B

S
M
L
XL
XXL

White

22341W
22343W
22345W
22347W

S
M
L
XL

Natural

22341N
22343N
22345N

S
M
L

Pink Fusion


22341PF
22343PF
22345PF

S
M
L

Reflector Yellow

22341RY
22343RY
22345RY
22347RY

S
M
L
XL





RELIEVES PATELLA TENDONITIS
EASES KNEE PAIN
REDUCES SWELLING
IMPROVES BALANCE

TARGETED CONDITIONS

Arthritis
Runners/Walkers Knee
ITBS (Iliotibial Band Syndrome)
Patellar Tendonitis
Patella Tracking

The Patented KS7 Knee Bracing Sleeve uses Compression Zone Technology® to embrace the knee joint in seven zones of graduated compression specifically designed to improve circulation while stabilizing the structure of the joint and relieving pressure on swollen and sore muscles and tendons. If you're facing runner's knee, jumper's knee, or the frustration of chronic knee pain, this thin and comfortable sleeve can offer the same benefits as bulkier braces, *and* you can wear it under your active wear or clothes.



ZONE 7
Light Compression
Wide, smooth & comfortable

ZONE 6
Moderate Compression
For increased muscle stability

ZONE 5
Firm Compression
"Knee Stabilizer Zone" for injury protection

ZONE 4
Light Compression
"Pressure Relief Zone"

ZONE 3
Firm Compression
"Lower Patella Stabilizer Zone" for pain relief & injury protection

ZONE 2
Moderate Compression
Boost venous flow and stabilize

ZONE 1
Light Compression
Wide, smooth, flat opening



Hypoallergenic gel
keeps sleeve from
slipping during activity



Where to measure	Size	Above Patella Circumference	Below Patella Circumference
	S	12-16 in/30-40 cm	9-13 in/23-33 cm
	M	14-18 in/35-45 cm	11-15 in/28-38 cm
	L	16-20 in/40-50 cm	13-17 in/33-44 cm
	XL	18.5-22.5 in/47-57 cm	15-20 in/38-50 cm
	XXL	20-25 in/53-63 cm	17-23 in/42-54 cm

Black

57340B
57342B
57344B
57346B
57348B

S
M
L
XL
XXL





ALL-IN-ONE QUAD AND ITB BRACE
RELIEVES ILIOTIBIAL BAND SYNDROME
EASES QUAD PAIN
REDUCES SWELLING

The Patent Pending QS4 Thigh Bracing Sleeve is the first, one-of-a-kind quad/thigh sleeve with an ITB brace already built in. It provides medical grade orthopedic support and graduated compression to ease the pain of ITB syndrome, quad strains, and hamstring weakness. Compression Zone Technology® makes the QS4 effective while still being easy to take on and off, light and comfortable to wear, even under normal active wear clothing. Unlike similar ITB straps that rely on velcro or other fasteners, the QS4 is safe, stays in place during any activity and does not shift placement during activity.

TARGETED CONDITIONS

ITBS (Iliotibial Band Syndrome)
Quad Strains
Hamstring Weakness
Leg Fatigue




Hamstring Stabilizer Zone



Built-in ITB Brace



Anti-migration gel system
Hundreds of hypoallergenic gel dots stabilize and keep sleeve from slipping



ZONE 4
Light Compression
Wide, smooth & comfortable with anti-migration gel system

ZONE 3
Moderate Compression
Muscle Stability

ZONE 2
Firm Compression
"ITB Stabilizer Zone"

ZONE 1
Light Compression
Non-binding flex band

Where to measure	Size	Mid-Thigh Circumference
	S	15-17 in / 31-43 cm
	M	17.5-20 in / 44.5-51 cm
	L	20.5-23.5 in / 52-60 cm
	XL	24-27 in / 60-68.5 cm

Black

77340B	S
77342B	M
77344B	L
77346B	XL





RELIEVES PATELLA TENDONITIS
CIRCULATION SAFE
EASES KNEE PAIN
IMPROVES PATELLAR TRACKING

TARGETED CONDITIONS

Arthritis
Runners Knee
Patellar Tendonitis
Knee Pain

TARGETED CONDITIONS

Lateral Epicondylitis
Medial Epicondylitis
General Elbow Pain



RELIEVES ELBOW PAIN
CIRCULATION SAFE
EASES MUSCLE SORENESS
WEAR WHILE ACTIVE

The Patented PS3 Patella Bracing Sleeve provides medical grade orthopedic support and graduated compression to ease the pain of patellar tendonitis (also known as jumper's knee). Compression Zone Technology® makes the PS3 effective while still being easy to take on and off, light and comfortable to wear – even under normal active wear clothing. Unlike similar patella straps that rely on velcro or other fasteners, the PS3 is safe, stays in place during any activity & will not bind behind the knee.



Hypoallergenic gel keeps sleeve from slipping during activity



ZONES

Where to measure	Size	Below Patella Circumference
	S	8-10.5 in/20-26 cm
	M	11-13.5 in/27-34 cm
	L	14-16.5 in/37-42 cm
	XL	17-20 in/43-50 cm



Where to measure	Size	Below Elbow Circumference
	S	8-10 in/20-25 cm
	M	10.5-12 in/26-31 cm
	L	12.5-15 in/32-38 cm
	XL	15.5-18 in/40-45 cm



The Patented ES3 Elbow Bracing Sleeve combines medical grade orthopedic support with the circulatory safety of graduated compression to treat painful lateral epicondylitis (also known as tennis elbow) and medial epicondylitis (or golfer's elbow). Using Compression Zone Technology®, the ES3 wraps the elbow joint in three graduated zones of compression to stabilize the tendons for relief of pain in any activity. Wear the ES3 under your normal shirt or jacket for all day relief and support.



ZONES



ZONE 3
Light Compression
Graduated compression for added comfort & increased flexibility

ZONE 2
Firm Compression
Firm Medical Grade compression to support & lift the patella

ZONE 1
Light Compression
Gentle compression with smooth, flat opening for non-binding fit


Black
72340B S
72342B M
72344B L
72346B XL

Black
62340B S
62342B M
62344B L
62346B XL

ZONE 3
Light Compression
Graduated compression for added comfort & increased flexibility

ZONE 2
Firm Compression
Firm Medical Grade compression to stabilize the damaged tendons

ZONE 1
Light Compression
Gentle compression with smooth, flat opening for non-binding fit





RELIEVES CARPAL TUNNEL
IMPROVES CIRCULATION
EASES WRIST & THUMB PAIN
REDUCES SWELLING

TARGETED CONDITIONS

Carpal Tunnel Syndrome
Arthritis
Wrist Pain
Wrist/Hand Fatigue
Swelling

The Patented WS6 Wrist Bracing Sleeve is a one-of-a-kind wrist brace designed to stabilize without immobilizing for relief while active. It provides medical grade orthopedic support and graduated compression to ease the pain of carpal tunnel syndrome, arthritis, sprains or general wrist pain. Hypoallergenic gel support system combined with Compression Zone Technology® keeps the sleeve in position while providing extra support and stabilization where needed most. Easily slips on and off. Thin profile fits under any apparel and sleeves. More comfortable and effective than traditional straps or cuffs. Designed to fit individual size. This is the most comfortable wrist support anywhere.



ZONE 6

Light Compression

With smooth, non-binding opening

ZONE 5

Moderate Compression

To relieve pressure and stimulate circulation

ZONE 3 & 4

Firm Compression

To support aching tendons & muscles, and further enhance circulation for pain relief

ZONE 2

Moderate Compression

For continued circulation support & reduces edema

ZONE 1

Light Compression

With smooth, non-binding opening



Hypoallergenic gel keeps sleeve from slipping during activity


Where to measure	Size	Wrist Circumference
	S	5-6 in/13-15 cm
	M	6.25-7 in/16-18 cm
	L	7.5-8.25 in/19-21 cm
	XL	8.5-9.25 in/22-24 cm

Black

82340B S
82342B M
82344B L
82346B XL

Natural

82340N S
82342N M
82344N L
82346N XL



The OrthoSleeve Retail Program

OrthoSleeve sets your business apart with a range of orthopedic supports that are one-of-a-kind and designed specifically for the needs & conditions of clients who wish to remain active & alleviate pain. OrthoSleeve helps you help your customers and become a comfort & wellness destination.

Retail Packaging

OrthoSleeve retail packaging was custom designed for the convenience of both the retailer and consumer. Bright colors & images help distinguish each OrthoSleeve product. Features and benefits are easy for staff and consumers to find and understand. A sticker at the top of the OrthoSleeve box on the package front allows the consumer to see the product and its color without opening the package. Indicated conditions are prominently displayed on the front to help determine which product is most suitable, and clear sizing instructions are on the side. Packages are sized appropriately for these lightweight bracing sleeves requiring less space for more profit.

Merchandising Flexibility

A variety of merchandising options are available to meet your space and offering requirements. OrthoSleeve offers easy start-up planograms of various sizes. Our packages are easily displayed on pegs and in a custom counter top/shelf display. A freestanding floor display rack is also available when wall space is at a premium. OrthoSleeve signage attracts customers to the profitable braces & supports section of your retail space.

Partnering Support

At OrthoSleeve, we understand the unique challenges faced by our retail partners. In addition to being recommended by Podiatrists nationwide and high levels of consumer satisfaction and loyalty amongst OrthoSleeve users, we work to grow brand recognition through advertising, active social media presence and event participation. OrthoSleeve retailers enjoy the benefit of being listed on our website dealer locator found at OrthoSleeve.com.

Store Profit

- Low space requirement
- High profit per square foot
- Patient compliance
- High levels of customer satisfaction
- Repeat sales & higher returns



Merchandising Aids



2' Slat Wall Display



Joint Pain Information Sign



Sell Sheets




Shelf Talkers



Slat Wall/Rack Header Sign



18"x24" Poster



Consumer Brochure



Prescription Pads



Tape Measure



Mannequin Foot



12 Pack P.O.P. Display



Freestanding Floor Display



Floor Display Art Panels

To Order Call 877-647-0386



OrthoSleeve® Family of Products



Compression Levels:

	ZONE 1	ZONE 2	ZONE 3	ZONE 4	ZONE 5	ZONE 6	ZONE 7	ZONE 8	ZONE 9	ZONE 10	ZONE 11
FS6	10-15mmHg	15-20mmHg	20-30mmHg	15-20mmHg	20-30mmHg	10-15mmHg	X	X	X	X	X
AF7	10-15mmHg	15-20mmHg	20-30mmHg	15-20mmHg	20-30mmHg	15-20mmHg	10-15mmHg	X	X	X	X
CS6	10-15mmHg	20-30mmHg	20-30mmHg	15-20mmHg	15-20mmHg	10-15mmHg	X	X	X	X	X
FS6+	10-15mmHg	15-20mmHg	20-30mmHg	15-20mmHg	20-30mmHg	10-15mmHg	20-30mmHg	15-20mmHg	10-15mmHg	15-20mmHg	10-15mmHg
KS7	10-15mmHg	15-20mmHg	20-30mmHg	10-15mmHg	20-30mmHg	15-20mmHg	10-15mmHg	X	X	X	X
PS3	10-15mmHg	20-30mmHg	10-15mmHg	X	X	X	X	X	X	X	X
QS4	10-15mmHg	20-30mmHg	15-20mmHg	10-15mmHg	X	X	X	X	X	X	X
ES3	10-15mmHg	20-30mmHg	10-15mmHg	X	X	X	X	X	X	X	X
WS6	10-15mmHg	15-20mmHg	20-30mmHg	20-30mmHg	15-20mmHg	10-15mmHg	X	X	X	X	X

Condition Chart:

CONDITIONS	FS6	AF7	CS6	FS6+	KS7	PS3	QS4	ES3	WS6
Achilles Tendonitis	✓	✓	✓	✓					
Ankle Pain/Strain	✓	✓							
Ankle Instability		✓							
Arch Pain	✓			✓					
Arthritis Pain	✓	✓	✓	✓	✓	✓		✓	✓
Calf Cramps			✓	✓					
Carpal Tunnel Syndrome									✓
Floating Patella					✓	✓			
Golfers Elbow								✓	
Hamstring Pain							✓		
Heel Pain	✓			✓					
Iliotibial Band Syndrome							✓		
Patellar Tendonitis					✓	✓			
Plantar Fasciitis	✓			✓					
Quad Strain							✓		
Severs Disease	✓			✓					
Shin Splints			✓	✓					
Swollen Feet	✓			✓					
Tennis Elbow								✓	



FAQ's:

Q. **How do I put the sleeves on?**

A. OrthoSleeves are shaped to fit *your* body. The products should be positioned on the body as shown on the front of our packages. Images illustrating proper sleeve positioning can also be found on our website at OrthoSleeve.com

Q. **How long will my OrthoSleeve last?**

A. All medical grade compression garments are designed to last six (6) months with normal wash & wear. Even if they look & feel normal, the LYCRA compression may begin to loosen after six months causing reduced efficacy.

Q. **How long can I wear my OrthoSleeves?**

A. You can wear the OrthoSleeve for up to 8 hours. Remove for one to two hours before wearing them again.

Q. **Are there any adverse reactions or warnings?**

A. If you are diabetic, consult your doctor before using; if you experience any discomfort or pain while using, remove product immediately; not intended for use on broken or injured skin.



GLOBAL DISTRIBUTION



ORDERING INFORMATION

ORDER ONLINE: Email orders to order@ingsource.com

FAX ORDERS: (877) 635-1521

CALL TO ORDER: (877) 647-0386

WWW.ORTHOSLEEVE.COM



Manufactured by: ING SOURCE, INC.
1720 TATE BLVD SE HICKORY, NC 28602
877-647-0386 INFO@ORTHOSLEEVE.COM

